

clinell[®]

Universal Wipes 100

Simple cleaning
& disinfection,
proven in practice

- ✓ Kills 99.99% of bacteria and enveloped viruses
- ✓ Effective from 30 seconds
- ✓ The UK's number 1 hospital wipe



Results you can trust

Where to use

Perfect for use on large surfaces and to tackle tough cleaning jobs. Suitable for offices, workspaces and high touch points.

Who are Clinell?

Since 2006, Clinell has been at the forefront of infection prevention. Today the Clinell Product Range has a wide selection of products to help clean and disinfect surfaces, but it all started with the Clinell Universal Wipe. Trusted by 9 out of 10 UK hospitals, Clinell Universal Wipes have become the go-to infection prevention solution within the healthcare industry.

Why does efficacy matter?

Many different disinfection products show how effective they are at removing microorganisms with a percentage on their packaging. This is called a log reduction. Log reductions are calculated by testing disinfectants against different microorganisms and seeing how well they kill them.

For example, if you put disinfectant in a petri dish with 1,000,000 microorganisms and it leaves 1,000 microorganisms, then that disinfectant would have a 3 log rating: equalling 99.9%. However, Clinell Universal Wipes has a 4 log rating, meaning that out of 1,000,000 microorganisms, only 100 remained, resulting in a 99.99% efficacy rating.

What about contact time?

Contact time determines how long a disinfectant must remain on a surface to reach its log reduction. Other disinfectants need 2-5 minutes and sometimes even 10 minutes to effectively disinfect a surface. Long contact times aren't practical as the surface must remain wet the entire time.

That's why fast contact times like that of Clinell Universal Wipes deliver an effective clean very quickly, making cleaning and disinfection simple.

The key to Clinell

Clinell Universal Wipes 100


Clinell Universal Wipes are designed for real-world applications, not just a petri dish. A lot of science goes into each wipe to provide fast and effective results you can trust, killing 99.99% of bacteria and enveloped viruses from 30 seconds.



Effective from 30 seconds



Kills 99.99% of bacteria and enveloped viruses




The UK's number 1 hospital wipe

How to use Clinell Universal Wipes




Risk assess the need for PPE.



Remove one wipe from the pack.



Working from top to bottom, clean to dirty, wipe in an 'S' shaped pattern.



Should your wipe dry out or become soiled, replace it with a new wipe. Let the surface air dry to maximise contact time.

Clinell Universal 100
is made by GAMA Healthcare.
For more information,
speak to your GAMA Healthcare
Area Manager or visit
www.gamahealthcare.com

Product	Unit of issue
Clinell Universal Wipes	Pack of 100 wipes
Order code	
BCW100	

Use biocides safely. Always read the label
and product information before use.



GAMA Healthcare Ltd., The Maylands Building,
Maylands Avenue, Hemel Hempstead, HP2 7TG, UK.
T: +44 (0)20 7993 0030 E: info@gamahealthcare.com
www.gamahealthcare.com

JBN220373